

## ► Thermo Scientific Nalgene Narrow-Mouth Packaging Bottles with Closures

natural high-density polyethylene with polypropylene closures, bulk pack



**Nalgene™ HDPE Narrow-Mouth Packaging Bottles** are ideal for specialty chemicals, diagnostics, biologicals, reagents, adhesives and veterinary pharmaceuticals.

Narrow-mouth HDPE packaging bottles are excellent for storage, shipping and packaging of liquid reagents.

### details

- Injection blow-molded for uniform wall thickness and strength
- Made from high-purity, food-grade resins
- Tough plastic construction
- Polypropylene closures included
- Leakproof!

**Recommended for:** Storage and shipment of liquid reagents.

**Ordering Information:** Bottles and closures are packed in separate polybags within the same carton.

**Compliance:** Manufactured in compliance with ISO 13485:2003 quality requirements. Bottles and Closures meet the requirements of 21CFR177.1520 for food and beverage use, USP Class VI.

### Nalgene Narrow-Mouth Packaging Bottles

Cat. No.	Capacity, mL (oz.)	Closure Size, mm	Neck I.D., mm	O.D., mm	Height with Closure, mm	No. per Pack	No. per Case
312089-0001	30 (1)	20	14	34	61	1000	1000
312089-0002	60 (2)	20	14	39	86	1000	1000
312089-0004	125 (4)	24	18	50	104	500	500
312089-0008	250 (8)	24	18	61	133	250	250
312089-0016	500 (16)	28	21	73	171	125	125
312089-0032	1000 (32)	38-430	27	92	217	50	50

## ▶ Thermo Scientific Nalgene Wide-Mouth Packaging Bottles with Closures

natural high-density polyethylene with polypropylene closure, bulk pack



**Nalgene™ HDPE Wide-Mouth Packaging Bottles with polypropylene closures are injection blow-molded for strength and uniform wall thickness.**

Durable, leakproof<sup>f</sup> chemical resistant bottles designed for shipping and storage.

### details

- Wide mouth for easy filling and dispensing of reagents
- Tough plastic construction
- Made from high-purity, food-grade resins
- Leakproof<sup>f</sup>

**Ordering Information:** Bottles and closures are packed in separate polybags within the same carton.

**Compliance:** Manufactured in compliance with ISO 13485:2003 quality requirements. Bottles and Closures meet the requirements of 21CFR177.1520 for food and beverage use, USP Class VI.

### Nalgene Wide-Mouth Packaging Bottles

Cat. No.	Capacity, mL (oz.)	Closure Size, mm	Neck I.D., mm	O.D., mm	Height with Closure, mm	No. per Pack	No. per Case
312189-0001	30 (1)	28	21	34	63	1000	1000
312189-0002	60 (2)	28	21	38	86	1000	1000
312189-0004	125 (4)	38	28	51	100	500	500
312189-0008	250 (8)	43	33	61	131	250	250
312189-0016	500 (16)	53	44	73	168	125	125
312189-0032	1000 (32)	63	52	91	199	50	50

## ▶ Thermo Scientific Nalgene Wide-Mouth Packaging Bottles with Closures

natural polypropylene copolymer with polypropylene closure, bulk pack



**Nalgene™ PPCO Wide-Mouth Packaging Bottles with polypropylene closures are ideal for diagnostics, biologicals, specialty chemicals, reagents, adhesives and veterinary pharmaceuticals.**

Durable, leakproof, autoclavable bottles are designed for shipping and storage.

### details

- Wide mouth allows easy filling and dispensing of powders
- Injection blow-molded for uniform wall thickness and strength
- Made from high-purity, food-grade resins; no plasticizers, fillers or extenders
- Polypropylene closures included
- Leakproof<sup>f</sup>

**Ordering Information:** Bottles and closures are packed in separate polybags within the same carton.

**Compliance:** Manufactured in compliance with ISO 13485:2003 quality requirements. Bottles and Closures meet the requirements of 21CFR177.1520 for food and beverage use, USP Class VI.

### Nalgene Wide-Mouth Packaging Bottles

Cat. No.	Capacity, mL (oz.)	Closure Size, mm	Neck I.D., mm	O.D., mm	Height with Closure, mm	No. per Pack	No. per Case
312187-0001	30 (1)	28	21	34	63	1000	1000
312187-0002	60 (2)	28	21	39	86	1000	1000
312187-0004	125 (4)	38	28	50	99	500	500
312187-0008	250 (8)	43	33	61	131	250	250
312187-0016	500 (16)	53	44	73	168	125	125
312187-0032	1000 (32)	63	53	91	199	50	50